APPLICATION OF SEMI-FITTED TOP SHEET


Note that there are two seams at each of the two bottom corners on the semi-fitted top sheet:

- The first seam is the corner seam
- The second seam (pleated seams) is where extra fabric has been added in for the foot comfort zone


INSTRUCTIONS:

1- Place the "corner seams" over the corners of the mattress

Note: We recommend that the seams be on the outside, so that they do not come into contact with the feet.



2- Pull the open end of the sheet up toward the head of the bed, stretching the material to its full length.



3- Return to the foot of the bed and find the "pleated seams". Grasp the seams and pull them toward the foot of the bed.

This will create a "pocket" at the foot of the bed that will allow for the feet of the patient/resident to be comfortably covered and free from pressure.

It is also helpful if you leave the blanket/warmspread untucked at the foot of the bed, as well.

