beurer

BM 54

CE 0483

DE	Blutdruckmessgerät Gebrauchsanweisung	2
EN	Blood pressure monitor Instructions for use	18
FR	Tensiomètre Mode d'emploi	32
ES	Tensiómetro Instrucciones de uso	48
IT	Misuratore di pressione Istruzioni per l'uso	64
TR	Tansiyon ölçme cihazı Kullanım kılavuzu	79
RU	Прибор для измерения кровяного давления в плечевой артерии Инструкция по применению	94
PL	Ciśnieniomierz Instrukcja obsługi	112
	Electromagnetic Compatibility Information	128

Inhalt

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Germany

Content

1. Getting to know your device	18
2. Important notes	19
3. Device description	21
4. Preparing the measurement	22
5. Measuring blood pressure	24
6. Evaluating results	26

7. Display and delete measured values	27
8. Transferring measured values	28
9. Cleaning and storing the device and cuff	29
10. Error messages/troubleshooting	29
11. Technical specifications	30
12. Warranty/service	31

Included in delivery

- Blood pressure monitor
- Upper arm cuff
- 4 x 1.5 V LR03 AAA batteries
- Storage bag
- Instructions for use

Dear Customer,

Thank you for choosing one of our products. Our brand stands for well-engineered, quality made products; driven by innovation and safety. Our product range encompasses a wealth of products for health and well-being in the heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air sectors. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the notes they contain.

With kind regards,
Your Beurer team

1. Getting to know your device

Check that the packaging of the Beurer BM 54 blood pressure monitor has not been tampered with and make sure that the required contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

The upper arm blood pressure monitor is used to carry out non-invasive measurement and monitoring of the arterial blood pressure values in adults.

This allows you to quickly and easily measure your blood pressure, save the measured values and display the development and average values of the measured values taken.

You are also warned of possible existing cardiac arrhythmia. The recorded values are classified and evaluated graphically.

2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Dispose of packaging in an environmentally friendly manner
	Manufacturer
	Permissible storage and transport temperature and humidity

	Permissible operating temperature and humidity
IP21	Protected against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling drops of water
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical devices.
	Certification symbol for products that are exported to the Russian Federation and members of the CIS

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measured values taken by you are for your information only – they are no substitute for a medical examination. Discuss the

measured values with your doctor and never base any medical decisions on them (e.g. medicines and their administration).

- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who is responsible for their safety or are instructed by such a person in how to use the device. Supervise children around the device to ensure they do not play with it.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the device on people who have the specified upper arm measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.

- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.
- Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty, the blood pressure monitor loses the date and time.
- To conserve the batteries, the blood pressure monitor switches off automatically if no buttons are pressed for 30 seconds.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.

- Only use the cuff included with the delivery or original replacement parts. Otherwise incorrect measured values will be recorded.

Notes on handling batteries


- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
-  **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush the batteries.

Instructions for repairs and disposal

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers.


- Do not open the device. Failure to comply will invalidate the warranty.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal. 

3. Device description

1. Cuff
2. Cuff line
3. Cuff connector
4. Display
5. Memory buttons **M1/M2**
6. START/STOP button **I**
7. Risk indicator
8. Connection for cuff connector (left-hand side)

Information on the display:

1. Time/date
2. Systolic pressure
3. Diastolic pressure
4. Calculated pulse value
5. Cardiac arrhythmia symbol
6. Pulse symbol
7. User memory
8. Risk indicator
9. Number of memory space/memory display for average value (P), morning (P^M), evening (P^M)
10. Release air (arrow)
11. Battery display symbol
12. Symbol for Bluetooth® transfer

System requirements for the “beurer HealthManager” app

- iOS ≥ 8.0 / Android™ ≥ 4.4
- Bluetooth® ≥ 4.0


List of compatible devices:

4. Preparing the measurement

Inserting the batteries


- Remove the battery compartment lid on the rear of the device.
- Insert four 1.5 V AAA (alkaline type LR03) batteries. Make sure that the batteries are inserted the correct way round. Do not use rechargeable batteries.
- Close the battery compartment lid again carefully.

All display elements are briefly displayed, 24 h or 12 h flashes in the display. Now set the date and time as described below. If the battery replacement symbol is permanently displayed, you can no longer perform any measurements and must replace all batteries. Once the batteries have been removed from the device, the date and time must be set again. Any saved measured values are retained.

Battery disposal

- The empty, completely discharged batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
- The codes below are printed on batteries containing harmful substances:
Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury

Set the hour format, date, time and Bluetooth® settings

The following section describes the functions and settings available on the blood pressure monitor.

Hour format → **Date** → **Time** → **Bluetooth®**

It is essential to set the correct date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

-  If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

Press and hold the START/STOP button for 5 seconds.

Hour format

- The hour format now flashes on the display.
- Select the desired hour format using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

Date

- The year flashes on the display.
- Select the year using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The month flashes on the display.

- Select the month using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The day flashes on the display.

- Select the current day using the **M1/M2** memory buttons and confirm with the **START/STOP** button .
-  If you have set the 12h hour format, the month is displayed before the day.

Time

The hours flash on the display.

- Select the current hours using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The minutes flash on the display.

- Select the current minutes using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

Bluetooth®

The **Bluetooth®** symbol flashes on the display.

- Use the **M1/M2** memory buttons to select whether automatic **Bluetooth®** data transfer is activated (**Bluetooth®** symbol flashes) or deactivated (**Bluetooth®** symbol is not shown) and confirm with the **START/STOP** button .


-  **Bluetooth®** transfers will reduce the battery life.

5. Measuring blood pressure


Ensure the device is at room temperature before measuring. The measurement can be performed on the left or right arm.

Attaching the cuff


Place the cuff onto the bare upper arm. The circulation of the arm must not be hindered by tight clothing or similar.


The cuff must be placed on the upper arm so that the bottom edge is positioned 2 – 3 cm above the elbow and over the artery. The line should point to the centre of the palm.


Now tighten the free end of the cuff, but make sure that it is not too tight around the arm and close the hook-and-loop fastener. The cuff should be fastened so that two fingers fit under the cuff.

Now insert the cuff line into the connection for the cuff connector.

This cuff is suitable for you if the index mark (▼) is within the OK range after fitting the cuff on the upper arm.


- ❗ If the measurement is performed on the right upper arm, the line should be located on the inside of your elbow. Ensure that your arm is not pressing on the line.

Blood pressure may vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm.

If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.

Important: The unit may only be operated with the original cuff. The cuff is suitable for an arm circumference of 22 to 44 cm.

Adopting the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying down. Always make sure that the cuff is at heart level.
- To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button ① to start the blood pressure monitor. All display elements are briefly displayed.

The blood pressure monitor will begin the measurement automatically after 3 seconds.

- ⓘ You can cancel the measurement at any time by pressing the **START/STOP** button ①.

As soon as a pulse is found, the pulse symbol ♥ will be displayed.

- The systolic pressure, diastolic pressure and pulse rate measurements are displayed.

- Er appears if the measurement could not be performed properly (see chapter 10 “Error messages/troubleshooting”). Repeat the measurement.

- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant M1 or M2 symbol appears on the display.

- Press the **START/STOP** button ① to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.
- If the device is not switched off manually, it will switch off automatically after 30 seconds.

If **Bluetooth® data transfer** is activated, data is transferred to the “beurer HealthManager” app after confirming the user memory with the **START/STOP** button ① (see chapter 8 “Transferring measured values”).

- The **Bluetooth®** symbol on the display flashes and the blue LED lights up. The blood pressure monitor now attempts to connect to the app for approx. 30 seconds.
- The **Bluetooth®** symbol stops flashing as soon as a connection is established. All measurement data is automatically transferred to the app. Once the data has been successfully transferred, the device switches off. If the data transfer was unsuccessful, the blue LED goes out and “Er” appears on the display.
- If a connection to the app cannot be established after 30 seconds, the **Bluetooth®** symbol goes out and the blood pressure monitor switches off automatically after 3 minutes.

- ⓘ Please note that you must add the blood pressure monitor to “My devices” in the “beurer HealthManager” app to enable data transfer. The “beurer HealthManager” app must be activated to allow data transfer. If the latest data is not displayed on your smartphone, repeat the data transfer as described in chapter 8.

If you forget to turn off the blood pressure monitor, it will switch off automatically after approx. 30 seconds. In this case too, the value is stored in the selected or most recently used user memory and the data is transferred if *Bluetooth®* data transfer is activated.

- Wait at least 5 minutes before taking another measurement.

6. Evaluating results

Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is a condition in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep. Arrhythmia can only be determined through an examination by your doctor. If the symbol  is shown on the display after the measurement has been taken, repeat the measurement. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor.

Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

Risk indicator:

The measurements can be classified and evaluated in accordance with the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Setting 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Setting 2: moderate hypertension	160–179	100–109	seek medical attention
Setting 1: mild hypertension	140–159	90–99	regular monitoring by doctor

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	<120	<80	self-monitoring

Source: WHO, 1999 (World Health Organization)

7. Display and delete measured values

User memory

The results of every successful measurement are stored together with the date and time. The oldest measurement is overwritten in the event of more than 60 measurements.

- Select the desired user memory ( ) with the **M1** or **M2** memory buttons if the device is switched off.
- To view the measurement data for user memory , press the **M1** memory button.
- To view the measurement data for user memory , press the **M2** memory button.
- The average of all measurements appears on the display.

User memory

- If *Bluetooth®* is activated (the  symbol flashes on the display), the blood pressure monitor attempts to connect to the app. If you press the **M1** button in the meantime, the transfer is cancelled and the average values are displayed. If you press the **M2** button, the transfer is cancelled and the measurement data of user memory  is displayed. The  symbol is no longer shown. The buttons are deactivated as soon as a connection is established and the data is transferred.

- ① Press the **M1** button if you have selected user memory 1.
Press the **M2** memory button if you have selected user memory 2.

A flashes on the display.

The average value of all saved measured values in this user memory is displayed.

- Press the relevant memory button (**M1** or **M2**).

AM flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

- Press the relevant memory button (**M1** or **M2**).

PM flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

Average values

- When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).
 - When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.
 - To switch the device off again, press the **START/STOP** button ①.
- ① You can exit the menu at any time by pressing the **START/STOP** button ①.

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start the retrieval of the average measured values. $\bar{x}$ flashes on the display and the average value of all saved measured values in this user memory is displayed.
- Press and hold the memory button **M1** or **M2** for 5 seconds depending on the user memory you are in.

All the values in the current user memory are deleted.

- To clear individual measurements from the relevant user memory, you must first select a user memory.
- Start the retrieval of the individual measured values.
- Press and hold the memory button **M1** or **M2** for 5 seconds (depending on the user memory you are in).
- The selected value is deleted. The device briefly displays CL 00.
- If you would like to delete other values, repeat the process described above.

You can switch off the device at any time by pressing the **START/STOP** ① button.

8. Transferring measured values

Transfer via *Bluetooth*® low energy technology

It is also possible to transfer the measured values saved on the device to your smartphone using *Bluetooth*® low energy technology.

You will need the “beurer HealthManager” app for this. This is available in the Apple App Store and from Google Play.

Proceed as follows to transfer values:

If *Bluetooth*® is activated in the settings menu, the data is transferred automatically after the measurement. The $\mathcal{B}$ symbol appears in the top left of the display (see chapter 4 “Preparing the measurement”).

Step 1: BM 54

Activate **Bluetooth®** on your device (see chapter 4 “Preparing the measurement”, **Bluetooth®**).

Step 2: “beurer HealthManager” app

In the “beurer HealthManager” app, add the BM 54 under “Settings/My devices”.

Step 3: BM 54

take a measurement.

Step 4: BM 54

data transfer immediately following measurements.

- If **Bluetooth®** data transfer has been activated, data is by transferred after having confirmed the user memory by pressing the pressing the **START/STOP** button ①.

Step 4: BM 54

Data transfer at a later point:

- Go to memory mode (chapter 7). Select the desired user memory. The **Bluetooth®** transfer starts automatically.

① The “beurer HealthManager” app must be active to allow data transfer.

If your smartphone has a protective cover, remove this to ensure that there is no interference during the transfer. Begin the data transfer in the “beurer HealthManager” app.

9. Cleaning and storing the device and cuff

- Clean the device and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances hold the device and cuff under water, as this can cause liquid to enter and damage the device and cuff.
- If you store the device and cuff, do not place heavy objects on the device and cuff. Remove the batteries. The cuff line should not be bent sharply.

10. Error messages/troubleshooting

In the event of errors, the error message $E_{r_}$ appears on the display.

Error messages may appear if:

- it was not possible to record the pulse: E_{r1} ;
- you move or speak during the measurement: E_{r2} ;
- the cuff is fastened too tightly or loosely: E_{r3} ;
- errors occur during the measurement: E_{r4} ;
- the pump pressure is higher than 300 mmHg: E_{r5} ;
- the batteries are almost empty  : E_{r6} ;
- E_{r7} the data could not be sent via **Bluetooth®**.

In such cases, repeat the measurement and/or data transfer. Ensure that the cuff line is correctly attached and that you do not move or speak during the measurement.

11. Technical specifications

Model No.	BM 54
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range	Cuff pressure 0 – 250 mmHg, Systolic 50 – 250 mmHg, Diastolic 30 – 200 mmHg, Pulse 40– 180 beats/minute
Display accuracy	Systolic ± 3 mmHg, Diastolic ± 3 mmHg, Pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: Systolic 8 mmHg / Diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 139 mm x W 94 mm x H 48 mm
Weight	349 g (without batteries, with cuff)
Cuff size	22 to 44 cm
Permissible operating conditions	+5°C to +40°C, 15-93% relative humidity (non-condensing)
Permissible storage conditions	-25°C to +70°C, $\leq 93\%$ relative humidity, 700– 1060 hPa ambient pressure
Power supply	4x 1.5 V --- AAA batteries

Battery life	Lasts for approx. 200 measurements depending on the blood pressure and pump pressure as well as the number of <i>Bluetooth</i> ® connections.
Accessories	Cuff, instructions for use, storage pouch
Classification	Internal supply, IP21, no AP or APG, continuous operation, application part type BF
Data transfer via <i>Bluetooth</i> ® wireless technology	The blood pressure monitor uses <i>Bluetooth</i> ® low energy technology, 2402MHz – 2480MHz frequency band, Transmission power max. 20 dBm, Compatible with Bluetooth 4.0 smart-phones/tablets

The serial number is located on the device or in the battery compartment.

Technical information is subject to change without notification to allow for updates.

- This device complies with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Services address or found at the end of the instructions for use.
- The device complies with the EU Medical Devices Directive 93/42/EEC, the German Medical Devices Act (Medizinproduktegesetz) and the standards EN1060-1 (non-invasive sphygmomanometers – Part 1: General requirements), EN1060-3

(non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC80601-2-30 (Medical electrical equipment – Part 2 – 30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).

- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.
- We hereby guarantee that this product complies with the European RED Directive 2014/53/EU. Please contact the specified service address to obtain further information, such as the CE Declaration of Conformity.
- Please contact the specified service address to obtain more detailed information such as the CE conformity declaration.

12. Warranty/service

In case of a claim under the warranty please contact your local dealer or the local representation which is mentioned in the list “service international”.

In case of returning the unit please add a copy of your receipt and a short report of the defect.

The following warranty terms shall apply:

1. The warranty period for BEURER products is either 5 years or- if longer- the country specific warranty period from date of Purchase.

In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.

2. Repairs (complete unit or parts of the unit) do not extend the warranty period.
3. The warranty shall not be valid for damages because of
 - a. improper treatment, e.g. nonobservance of the user instructions.
 - b. repairs or tampering by the customer or unauthorised third parties.
 - c. transport from the manufacturer to the consumer or during transport to the service centre.
 - d. The warranty shall not be valid for accessories which are subject to normal wear and tear (cuff, batteries etc.).
4. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

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Electromagnetic compatibility

Electromagnetic emissions

Guidance and manufacturer's declaration – electromagnetic emissions for all EQUIPMENT and SYSTEMS

The Blood pressure monitor (BM 54) is intended for use in the electromagnetic environment specified below. The customer of the user of the Blood pressure monitor (BM 54) should assure that it is used in such and environment.

Emission test	Compliance	Electromagnetic environment–guidance
RF emissions CISPR 11	Group 1	The Blood pressure monitor (BM 54) uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The Blood pressure monitor (BM 54) is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Electromagnetic immunity

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS

The Blood pressure monitor (BM 54) is intended for use in the electromagnetic environment specified below. The customer of the user of the Blood pressure monitor (BM 54) should assure that it is used in such and environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment–guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

The Blood pressure monitor (BM 54) is intended for use in the electromagnetic environment specified below. The customer or the user of Blood pressure monitor (BM 54) should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Blood pressure monitor (BM 54), including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance.
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz

			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Blood pressure monitor (BM 54) is used exceeds the applicable RF compliance level above, the Blood pressure monitor (BM 54) should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Blood pressure monitor (BM 54). ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT or SYSTEM that are not LIFESUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the Blood pressure monitor (BM 54).

The Blood pressure monitor (BM 54) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Blood pressure monitor (BM 54) can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Blood pressure monitor (BM 54) as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.1167	0.2334
0,1	0.3689	0.7378
1	1.1667	2.3334
10	3.6893	7.3786
100	11.6667	23.3334

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Android is a trademark of Google LLC.

Beurer GmbH • Söflinger Straße 218 • 89077 Ulm, Germany • www.beurer.com
www.beurer-blutdruck.de • www.beurer-gesundheitsratgeber.com • www.beurer-healthguide.com

655.12_BM54_2018-02-27_01_IM1_BEU