Please do not operate or allow anyone else to operate the Device without reading this User Manual.

Please use the Device only as directed.

www.arcmicrotech.com
www.arc4health.com

ARC Microtech Limited

Unit 7, Woodland Enterprise Centre, Hastings Road
Flimwell, East Sussex TN5 7PR, UK

© 2021 ARC Microtech Limited
Patent Numbers GB 2502133 & GB 2534681
Ver. 8.5

User Manual

The complete pain management
and tissue repair system

What is the Arc4Health Device?

Arc4Health is a remarkably small Device that is wearable, drug free and non-invasive. It works by automatically generating a complex series of microcurrents with a simple, one finger operation. In fact, it delivers over 4.6 million permutations within a single programme to support the body's own natural repair process. And just like naturally occurring microcurrents, they can't be felt, so cause no discomfort to you. Arc4Health has regulatory approval as a Class IIa Medical Device, meaning it is safe to use on humans. It is supplied subject to our terms and conditions of purchase available on our website. Please ensure that you read these prior to using the Device.

What's in the box?

- The Device • Arm Strap • Conductive Pads (x4)
- Carbon Rubber Pads (x2) • Conductive Gel (250ml)
- Charger • User Manual

Answers to most frequently asked questions

- No, the Device cannot be placed over an open wound or skin irritation.
- Yes, the Arc4Health is safe to use on injuries that have not broken the skin, providing the wearer is comfortable doing this.
- Always rotate the Device between arms and legs for each 3-hour treatment.

- Like any electrical device, do not use in wet or damp conditions.
- Please do not operate or allow anyone else to operate the Device without reading this User Manual. Please use the Device only as directed.
- Intended user: general public, no special skills, training or prior knowledge required. It is expected that the patient will operate the device for use on themselves.

Contraindications

- Do not use the Device if you have a latex or rubber allergy as some of the component parts and consumables may contain latex or rubber.
- Do not use the Device on babies, infants or children.
- The safety of Microcurrent Therapy during pregnancy has not been established, therefore pregnant users shall follow precautions recommended by their physicians.
- The safety of Microcurrent Therapy with cancer has not been established. Users with suspected or diagnosed cancer of any type shall follow precautions recommended by their physicians.

- Users with suspected or diagnosed epilepsy, or heart rhythm disorder, shall follow precautions recommended by physicians.
- Users with suspected or diagnosed heart disease and heart conditions of any type shall follow precautions recommended by their physicians.
- Do not use the Device at the same time as a TENS machine.

Programmes

We always recommend you run through the following 6-week treatment protocol first in order to maximise results. Once the initial 6-week protocol finishes you can select the individual programme that best suits your condition (please see below for a description of each programme). Every few months, we recommend running through the 6-week protocol again in order to maintain optimum results (please note that this is a manufacturers guidance and not a clinical recommendation).

As a general rule¹: (i) acute phases of pain and tissue repair are defined as lasting less than 12 weeks; and (ii) chronic phases of pain and tissue repair are defined as lasting more than 12 weeks.

- Indicator lights
- ⓘ On/Off and Programme button
- ▶ Play button

¹ www.britishpainsociety.org/media-resources

6 Week Protocol:

P1 = Pain Management (Acute Phase). 1 programme daily for **2 weeks**.

P2 = Pain Management (Chronic Phase). 1 programme daily for **1 week**.

P3 = Tissue Repair (Acute Phase). 1 programme daily for **2 weeks**.

P4 = Tissue Repair (Chronic Phase). 1 programme daily for **1 week**.

Each programme lasts for 3-hours.
P = Programme

How to put the Device on


Please inspect the Device before use and do not use if the case is damaged or cracked.

Only apply the Device to unbroken, clean, healthy skin.


Front (device in pocket)


Back (sticky conductive pads)


Back (carbon rubber pads)

Conductive Pads: Fix to the arm strap using the press studs. Please ensure that the conductive pads are placed back onto the Perspex film after use to help prolong their adhesiveness.

Carbon Rubber Pads: Fix to the arm strap using the press studs. Carbon rubber pads should be wiped clean with warm water after each use.

Please note: the conductive pads and carbon rubber pads work in the same way, fitting onto the arm strap using the press studs. The conductive pads have a sticky conductive surface, whereas the carbon rubber pads will need to be used with the conductive gel. You decide which pads you prefer to use.

Conductive Gel: Apply a small amount to the carbon rubber pads. The gel ensures a good electrical connection with the skin and is Propylene Glycol free. The Conductive Gel only needs to be used with the Carbon Rubber Pads.

Please Note: Do not use any other gel product in conjunction with the Device.


Arm Strap: Attach the arm strap with the Device to anywhere comfortable on the arm or leg. Rotate the strap between a different arm or leg after each 3-hour treatment. Please do not over tighten - the carbon/conductive pads are very sensitive and only light pressure is required to maintain connectivity. After putting the Device into the arm strap's pouch, please ensure that the electrical lead in the arm strap is inserted fully into the bottom of the Device.

- Never strap the Device over an open wound or skin irritation.
- The Arc4Health is safe to use on injuries that have not broken the skin, providing you feel comfortable with this.
- Always rotate the Device between different arms or legs for each 3-hour treatment.

Programming the Device


On/Off left button.
Press and hold for 3 seconds
when turning on or off.


All four indicator lights will very quickly show green, then yellow to indicate that your Device is functioning correctly.


After this, lights will flash green corresponding with the battery charge available: 1 light = 25%, 2 = 50%, 3 = 75%, 4 = 100%.

After 3 seconds all lights will go out. Immediately after this, a yellow light will indicate which programme was last used. 1st light = **P1**, 2nd = **P2**, 3rd = **P3**, 4th = **P4**.


Press the left (On/Off) button quickly to skip through and select the desired programme from **P1** to **P4**.


Press and release the right play button quickly to activate the selected programme. The unit will beep on activation of the selected programme.

Once activated, the right play button will turn either green (connected) or yellow (not connected).

After 30 seconds the light will go out to conserve battery, but the programme indicator will continue to flash throughout the delivery cycle.

If not connected, ensure you have used sufficient conductive gel on the carbon rubber pads and the electrical lead from the arm strap is still attached to the bottom of the Device.

Once the 3-hour delivery cycle is complete, the Device will turn itself off.

Cleaning your Device

If you intend to share your device with another, please ensure this is only after appropriate inspection and thorough cleaning with disinfectant wipes. Ensure they read this User Manual prior to use and use the Device only as directed.


Please also clean the Device as detailed above between uses. The following accessories should not be shared, with each individual using or purchasing their own: Carbon Rubber Pads, Conductive Pads and Arm Straps.

These should not be shared in order to prevent cross-infection, as they are skin-contacting accessories.

Charging

Please inspect the state of the plug charger (DC PSU) before use.

Device must be connected directly to mains power (230v, 50Hz) using the plug charger provided for charging to take place.

Only use the charger supplied, or an original replacement charger ordered directly from Device Manufacturer.

On receiving the Device, please charge fully before initial use. One 4 - 5 hour charge can give up to 30 hours of use, enough for 10 delivery cycles.

Insert the charger lead into the base of the Device. Plug into power socket. Two beeps indicate charging.

Battery charge is indicated by the corresponding number of flashing green and yellow 'indicator lights':

1 yellow = charging 0-25%


1 green + 1 yellow = charging 25-50%

2 green + 1 yellow = charging 50-75%

3 green + 1 yellow = charging 75-100%

4 green = fully charged

All indicator lights off = 100% charged

Please Note:

The Arc4Health delivery unit uses a lithium-ion battery. It needs to be charged on a reasonably regular basis in order to maintain optimal battery life and condition. Avoid long periods without using or recharging the device.

Turning the audible beep on and off

Press and hold both buttons at the same time until you hear three short beeps (approx. 3-5 seconds):

- High pitch – High pitch – High pitch = audible beep ON.
- High pitch – High pitch – Low pitch = audible beep OFF.

General Warnings

- Use this Device only with the leads, electrodes and accessories recommended by the manufacturer, and do not interconnect with any other equipment. Use of other parts or materials not supplied by the manufacturer may prove hazardous and will invalidate the warranty. For example, using an unapproved charger may produce an electrical fault and over-heating of the Device.
- Use of the Device is in addition to, and in no way a substitute for your normal, existing practice of seeking help, guidance and support from your GP, or other Clinical Practitioner.
- Use of Microcurrent Devices may cause some dehydration. Please ensure you drink a sufficient amount of water and stay hydrated when using the Device.
- Do not immerse any part of the unit in water.
- Do not use the Device on people with cardiac pacemakers, implanted defibrillators, or any other implanted active electronic devices, as these may cause electrical interference. Likewise, the Arc4Health Device may cause interference with implanted electrical devices. Please do not use these in conjunction in order to minimise interference.
- There is the potential for electromagnetic interference or other interference (e.g. magnets and radio frequency), between the Device and other devices in the home environment. This includes mobile phones, Wi-Fi device, electromagnets and Bluetooth devices. To minimise interference between devices, please keep the Device away (e.g. 1m) from other electrical devices/appliances.
- Keep the Device at least 3 feet / 1 meter away from excess heat (e.g. fireplaces or radiant heaters) and excess moisture (e.g. nebulisers and steam kettles), as these may damage the internal electronics and interfere with the electrical output.
- Keep the Device away from sources of radiation, such a microwaves and mobile phones, as this may interfere with the electrical output. Operation in close proximity (e.g. 1m) to a shortwave or microwave therapy device may produce instability in the Arc4Health output.

- Please do not use simultaneously with high frequency surgical equipment.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Microcurrents should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.
- In between uses, please store the Device and its accessories in the box to prevent a build-up of lint and dust which may damage the Device. This will also keep it out of direct sunlight, which may cause the Device to overheat.
- Keep the Device out of reach of children, pests and pets to prevent damage being caused e.g. from chewing/biting.

WARNING – STRANGULATION HAZARD

The charging cable may be a strangulation hazard for small children which could result in injury or death. Keep the charging cable out of the reach of small children (more than 3 feet / 1 meter away) and do not allow them to play with it.

Troubleshooting

- Right activation button showing yellow = loss of connectivity. Please refer to the Good/bad connectivity notes.
- Left On/Off button showing yellow = lack of charge. Device unit needs charging. Please refer to the Charging notes.
- Device keeps switching off = the right activation button was not pressed.
- Audible beep has gone = Device's default setting is to have audible sounds. Please refer to the notes for turning the audible beep on and off.

Important

- Do not place or strap the Device over an open wound or skin irritation.
- Do not apply any other creams or gels to the carbon rubber pads other than ARC approved conductive gel.
- The Device is in addition to, and in no way a substitute for, your existing medical treatment or any medical advice you have received from your healthcare professional.
- Always consult a medical expert/healthcare professional on any health conditions and on using the Device prior to use.
- Like any electrical device, do not use in wet or damp conditions.
- Some people, after prolonged use, may develop a sensitivity condition. This can manifest in a 'rub' type condition. We recommend that the arm or leg being used is alternated after every 3-hour cycle.

- If this occurs, users should stop using the Device and should contact the Device Manufacturers if they experience any adverse reactions. Contact details can be found on the back page.
- Please use the Device only as directed.
- Although all accessible parts are compliant with the biocompatibility requirements of EN 10993-5 and EN 10993-10, there is the potential that some users may experience an allergic reaction to the materials used. If this occurs, users should stop using the Device and should contact the Device Manufacturers if they experience any adverse reactions. Contact details can be found on the back page.

Good/bad connectivity

If the right play button turns green, the programme has started and you have good connectivity.

If the right activation button turns yellow, the programme has started but you DO NOT have connectivity.

If YELLOW, please check that:

- the electrical lead from the arm strap is still attached to the bottom of the Device;
- the carbon/conductive pads have not come loose or fallen off;
- the conductive pads have sufficient adhesiveness and have not dried out;
- the carbon rubber pads have sufficient gel applied.

Specifications

Arc4Health Delivery Unit

Battery	Li-ion Rechargeable, 3.7V nominal (fully traceable) – non removable.
Device	Rated voltage 12V
Pad connect	Via 2.1 x 5.5mm (centre positive) plug on unit
Maximum pad voltage	+/- 9V
Programme selection	Via push buttons
Modes	4 programmes Pulse duration: 0.11msec Pulse to pulse spacing: 0.86msec Pulse polarity: bipolar
Output channels	Single channel
Output voltage	0-9V (+/-), controlled to maintain the specified current Max = 4V for R = 10kohm (typical) Max = 0.2V for R = 500 ohms
Output current	0-400 uA
Tolerances	All output parameters less than +/- 10%

Physical dimensions	45mm (W) x 15mm (D) x 105mm (L)
Approximate weight	65g
Operating temperature	5°C to 40°C
Storage temperature	-25°C to 70°C
Please allow the Device 30mins to warm up from the minimum storage temperature and cool down from the maximum storage temperature before using at an ambient temperature (20 °C).	
Relative humidity	15 to 90%
Atmospheric pressure	700 mBar to 1060 mBar
IP rating	IP32 (ingress protection) 3 – protected from tools and wires greater than 2.5mm 2 – protected from water spray when tilted less than 15 degrees from vertical
Expected life of unit	3 years under normal use

Arc4Health System

MDD Classification	Class IIa (Rule 9)
Assessment Route	Annex V
Ingress Protection	IP32
Sterilisation	Not suitable to be sterilised

Arm Strap

MDD Classification	Class I (Rule 1)	
Assessment Route	Annex VII	
Type BF Applied Part		
Compliant with EN 10993-5 and EN 10993-10		
Material Composition: Arm strap outer fabric – OK fabric (80% Nylon, 20% Spandex Fibre) Arm strap inner fabric – Woven fabric (99% polyester, 1% conductive yarn)		
Expected Service Life – up to 4 months with moderate use		
Expected Shelf Life – up to 24 months		

Carbon Rubber Electrode

MDD Classification	Class I (Rule 1)	
Assessment Route	Annex VII	
Type BF Applied Part		
Compliant with EN 10993-5 and EN 10993-10		
Expected Service Life – up to 12 months with moderate use		
Expected Shelf Life – up to 24 months		
Degraded electrodes can cause connectivity issues, which may negatively impact the performance of the Device. Please ensure you replace your Carbon Rubber Electrodes after 12 months of use.		

Self-Adhesive Electrode

MDD Classification	Class I (Rule 1)	
Assessment Route	Annex VII	
Type BF Applied Part		
Compliant with EN 10993-5 and EN 10993-10		
Expected Service Life – up to 4 weeks with moderate use		
Expected Shelf Life – please refer to the 'use by' date printed on the packaging		
Degraded electrodes can cause connectivity issues, which may negatively impact the performance of the Device. Please ensure you replace your Self-Adhesive Electrodes after 4 weeks of use.		

Conductive Gel

MDD Classification	Class I (Rule 1)
Assessment Route	Annex VII
Compliant with EN 10993-5 and EN 10993-10	
Expected Service Life – up to 4 weeks with moderate use	
Expected Shelf Life – please refer to the 'use by' date printed on the bottle.	

DC Power Supply (charger)

Only use the following supplied charger when charging the Device, or an original replacement charger ordered directly from Device Manufacturer.	
Connection	Via 2.1 x 5.5mm (centre positive) plug on unit (pads removed whilst charging)
Make	Dee Van Enterprise Co, Switching Adaptor
Model	DSA-12PF09-12 bc 120y
REF	ARCCHARGE235001
Charger input	100-240 VAC, 50-60Hz, 0.5A
Charger output	12 VDC, 1A
Compliance	EN 55032:2015, EN 61000-3-2:2014, EN 61000-3-3:2013, EN 55035:2017

DC Power Supply Symbols

	Charger should be operated indoors only, where it is not exposed to the weather.
	The Charger has approval as Class II equipment.
Expected Service Life – up to 12 months with moderate use	
Expected Shelf Life – up to 24 months	

Warranty Our limited warranty covers defects in materials and workmanship or electrical fault under normal use for 2 years. This is valid from the date of delivery of the Device to the original purchaser. If a defect arises and a valid claim is received within the warranty period, we will either refund the cost of the defective Device or exchange the Device with a Device that is new, or which has been manufactured from new or serviceable used parts.

WARNING: No modification of this equipment is allowed.

Please note: There are no user serviceable parts, the unit must be sent back to the manufacturer if a repair is required.

Exclusions: Carbon Rubber Pads, Conductive Pads, Conductive Gel, Arm Straps and Chargers are not covered by this warranty - **only the Device.**

This warranty is set out in and subject to our terms and conditions of purchase.

The Device complies with the requirements of EN 60601-1-2, but it is recommended that close proximity to portable and mobile RF communications such devices (Wi-Fi, hubs, mobile and cordless telephones etc) be restricted during treatments as these can influence medical electrical equipment.

WEEE Statement

The Symbol on the product indicates that it should not be disposed of with your normal household waste but should be taken to a designated collection point for the recycling of waste electrical and electronic equipment.

RoHS statement

The Device is subject to, and compliant with, the requirements of the Restriction on Hazardous Substances (RoHS) Directive, 2011/65/EU.

Standard Symbols

CE Marking

Manufacturer

Refer to instruction booklet

Type BF Applied Part

Single Patient Use

EU Declaration of Conformity

The Device Microcurrent Therapy System complies with the requirements of the Medical Devices Directive 93/42/EEC as amended. The Device is considered to be a Class IIa Medical Device.

Compliance has been demonstrated with reference to the following harmonised standards:

- EN 60601-1
- EN 60601-2-10
- EN 60601-1-11

Contact details

For further information please direct any general enquiries to:

admin@arcmicrotech.com

Tel: +44 (0) 1580 755504

Please note: All enquiries will be handled between 09:00 – 17:00 Monday - Thursday, 09:00 – 16:00 Friday.

Please visit **www.arc4health.com** or **www.arcmicrotech.com** should you need more information on our Device (including our FAQs) or if you need to buy more consumables. Please use the following details to contact us for assistance if needed, in setting up, using or maintaining the Device. Please report any unexpected operation or events to us immediately.

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