

# Product Data Sheet

2/4/2018

## FEATURES

- Three section swing arm offering height, rotation & angle adjustment
- Compatible with the PT range of single & double door pharmacy trolleys
- Laptop holder fitted to an upright bolted to the rear left hand corner of the trolley (upright & holder supplied for customer to fit)
- Adjustable to accept laptops 280-435 x 22mm (w x h)
- Laptop supported in a cradle fitted with a steel mounting bracket that can be fixed using either nylock or wing nuts
- Laptop can be secured using a security clamp (supplied)
- Design allows easy access to laptop ports
- Retrofit model available, compatible with existing PT pharmacy trolleys

## STANDARD COLOURS

- Paint - Grey

## PRODUCT LINK

[http://www.bristolmaid.com/product.asp?CAT\\_ID=0&P\\_ID=1008](http://www.bristolmaid.com/product.asp?CAT_ID=0&P_ID=1008)



Item	Description	Weight
PT/LHNR/S	Non-Removable, Swing Arm	5.0kg
R/PT/LHNR/S	Non-Removable, Swing Arm (Retro-Fit)	5.0kg

For information on recycling & cleaning instructions please refer to [www.bristolmaid.com](http://www.bristolmaid.com) or contact our technical services department on 01258 451338

Blandford Heights, Blandford Forum, Dorset DT11 7TG, UK

Tel: 01258 451338

Fax: 01258 455056

Email: [sales@bristolmaid.com](mailto:sales@bristolmaid.com)

[www.bristolmaid.com](http://www.bristolmaid.com)



**Bristol Maid**<sup>TM</sup>  
Hospital Metalcraft Limited 