

EC Declaration of Product Conformity

The manufacturer established in the community;

**BM Polyco Ltd
Crown Road
Enfield
EN1 1TX
United Kingdom**

Certified to ISO 9001 & ISO 13485

declares that the new Medical Device described hereafter



FINITY

Stretch Vinyl Powder Free Synthetic Disposable Glove

is in conformity with the provisions of Council Directive 93/42/EEC and with National Standards transposing harmonised standards EN455-1, EN455-2, EN455-3 and EN455-4

is Class I and is self-certified according to the requirements of EC Directive 93/42/EEC and subsequent amendments

Product Information

Glove Care: Store below 40°C away from direct sunlight

Freedom from holes: AQL 1.5

Physical properties: 3.6N minimum



User Information

This product does not contain natural rubber latex or chemical accelerators.

Range Available: sizes and codes

Sizes	5½	6½	7½	8½	9½
24cm	-	FT100/01	FT100/02	FT100/03	FT100/04
30cm	FT130/00	FT130/01	FT130/02	FT130/03	FT130/04

Done at Enfield, 16/01/14



Bernard Garvey
Technology Director