

**Get a good grip on handling
even the toughest suction scenario.**



V-VAC™ Manual Suction Unit

Get a good grip on handling even the toughest suction scenario with Laerdal's V-VAC Manual Suction Unit. It's ready to go, easy to use, and up to the task. Ideal for First Responder emergency kits and as back-up suction, it meets all your suction needs, infant through adult. One-handed operation, one-piece disposable cartridge, and unique no-clog suction tip make the V-VAC the manual suction unit of choice!



Laerdal®
helping save lives

Product Specifications

FEATURES	BENEFITS
Performance	<ul style="list-style-type: none"> All in one system adjusts for high and low vacuum suctioning Adapter tip allows for wide range of catheter connections
Efficiency	<ul style="list-style-type: none"> Exceeds 380 mmHg suction with peak air flow of over 70 liters per minute Unique no-clog suction tip resists clogging
Convenience	<ul style="list-style-type: none"> Disposable one piece collection cartridge makes cleanup easy Light weight design and compact size allow for storage in emergency bags, nursing offices, and first aid stations
Economy	<ul style="list-style-type: none"> Reusable handle is easy to clean and put back into repeated service Low cost and manual operation make V-VAC™ Manual Suction Unit the right choice for extended storage in mass casualty kits



The V-VAC unit fits easily into many emergency kits. A separate carry bag is also available.



The soft tip of the intake valve has deep cutouts to prevent tube occlusion and soft tissue injury. The tri-cut valve resists back flow of suctioned matter.

Part No.	Description
98 50 00	V-VAC Manual Suction Unit Starter Kit*
98 51 00	V-VAC Manual Suction Unit Training Kit*
98 50 01	Replacement Cartridge
98 50 30	V-VAC Handle
98 52 00	Carry Bag

*Available in cases of 15. Call for details.

V-VAC Starter Kit for Patient Use Includes:

Handle, 2 Replacement Cartridges, Short Suction Catheter w/Adapter Tip, Double Male Connector, and Directions for Use.

V-VAC Training Kit for Practical Use Includes:

Handle, 1 Training Cartridge with plug for emptying and cleaning, Catheter w/Adapter Tip, and Double Male Connector.



For additional information or to order, contact your local Authorized Laerdal Distributor!

Distributed in USA by:
 Laerdal Medical Corporation
 167 Myers Corners Road,
 Wappingers Falls, NY 12590
 (877) LAERDAL (523-7325)
 or (845) 297-7770
 FAX (800) 227-1143
 customerservice@Laerdal.com
 www.Laerdal.com

Distributed in Canada by:
 Laerdal Medical Canada Ltd.
 151 Nashdene Rd., Unit #45
 Toronto, ON M1V 4C3, CANADA
 (416) 298-9600 or (888) LAERDAL
 FAX (416) 298-8016
 savelives@Laerdal.ca