



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical sheet

Anti-Fog Solution with Sponge

Reference(s)

NONFB100

Description

Anti-Fog solution with Sponge for clear view
Sponge with adhesive on the back for easy application

Main materials

Solution: Dilute isopropyl alcohol 6ml, concentration 10%
Sponge: Polyester
Bottle: Polyethylene

Pack factor

20 units per carton

Packaging

Peel pouch

Manufacturer

Aspen Surgical Products, Inc.
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: IIa sterile