




Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant
www.medline.com/en/



Invisishield™

Reference(s)

IDE1015 – IDE1019

Description

References	A- Width (cm)	B- Length (cm)	Adhesive Split (cm)
IDE1015	122	132	15x48 cm
IDE1019	122	180	15x48 cm

Main materials

Clear Frosted PE film, 50 µm

Pack factor

40 units per carton
10 units per box – 4 boxes

Packaging

Paper and film peel pouch

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile



Latex free



0086



STERILE EO

© 2016 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_InvisishieldDrape_IDE1015-IDE1019_EN01

01.01.INT.12