



About Talar Made

TalarMade Ltd are based in the UK and employ a large number of healthcare professionals that form part of our management and clinical team. We pride ourselves on our ability to transfer our medical expertise and market knowledge into high quality products, which ultimately leads to improvements in patient care.

Following successful consultation with a wide range of UK leaders in the medical field of pressure redistribution, diabetic foot care, acute and community healthcare providers and tissue viability specialists we have developed a range of high quality devices that offer the clinician a number of solutions that may be incorporated into a pressure damage prevention strategy.

Our Pressure Care range includes:

•	Safeguard Heel Boot	Page 3
•	Safeguard Utility Pad	Page 3
•	Heelpro	Page 4
•	Heelpro Advance	Page 4
•	Footsafe Prevention Boot	Page 5
•	HeelSafe Pressure Redistribution Pad	Page 6
•	SoleSafe Bed End Pressure Redistribution Pad	Page 6
•	FootShield II Shoe	Page 7
•	FootShield II Boot	Page 7
•	FootShield II Wedge Shoe	Page 8
•	FootShield Hex Insole	Page 8
•	Ambulatory Pressure Relief Boot	Page 9
•	Pressure Relief AFO	Page 9

Our helpful product selection pathways are on pages 10-11 if you need assistance selecting any products within our range.

To order our products or for more information, you may contact us at the details below.



01246 268456



info@talarmade.com



talarmade.com



Features a wide, soft and secure fixation strap that is easy to place and remove.


Safeguard Heel Boot

The Safeguard Heel Boot has been designed to redistribute pressure on the bony prominences of the ankle, foot and heel.

It is designed to assist in the prevention of pressure ulcers as part of a patient's overall pressure ulcer prevention care program.

- * Supplied in pairs
- * Single patient use

Order Code	PIP Code	NHSSC Code	Qty
SGFB/U	413-2536	FES11322	1
SGFB/U10		FES11255	10



L 50 x W 30 x H 2.5 cm Product can be cut to shape

Safeguard

Utility Pad

The Safeguard Utility Pad is designed to assist in the prevention of pressure ulcers as part of a patient's overall pressure ulcer prevention care program.

It's designed redistributes pressure on bony prominences such as knees, elbows, spinous processes and heels.

- * Supplied in pairs
- * Single patient use

Order Code	PIP Code	NHSSC Code	Qty
SGFUP/U	413-2544	FES11251	1
SGFU/U10		FES11308	10



heelpro

For the prevention and management of heel pressure ulcers

- Completely offloads by elevating the heel
- Helps reduce incidents of heel pressure ulcers¹
- Heel aperture allows for easy inspection
- Universal sizing for easy application
- Vacuum packed to reduce storage space
- Available to prescribe on FP10
- * Supplied individually or in packs of 10
- * Single patient use

Order Code	NHSSC Code	PIP Code	Qty
HPRO/U	GRD1502	409-1476	1
HPRO/U10	GAA123		10



Generous low friction cushioning for effective protection

Replacement Lycra Panels

Order Code	NHSSC Code	Qty
HPRO/RLP	GRA1663	10



heelpro Stability Wedge

Order Code	NHSSC Code	Qty
HPRO/SW/U	GRA5182	1
HPRO/SW/U10	GRA5183	10

heelpro Advance

- Completely offloads by elevating the heel
- Wipe clean cover which means it is ideal for use where there are existing wounds
- Heel aperture allows for easy inspection
- Universal sizing for easy application
- Vacuum packed to reduce storage space
- * Supplied individually or in packs of 10
- * Single patient use

Order Code	NHSSC Code	PIP Code	Qty
HPRO-AD/U	GAA092	414-6411	1
HPRO/AD/U10	GAA124		10



FootSafe

Prevention Boot

- The FootSafe Boot may contribute to interventions in preventing and managing heel ulcers².
- Completely off-loads by elevating the heel
- Supplied inflatable or pre-inflated
- Anti-microbial material
- Easy close PU buttons and straps
- Durable

* Supplied individually or in packs of 10

* Multi patient use



Available in 5 sizes.

Style	Size	Qty	Product Code	NHSSC Code
Pre-Inflated	Child	1	PHPB-UC-PI-K-P	
Pre-Inflated	Child	10	PHPB-UC-PI-K-P5	
Pre-Inflated	Youth	1	PHPB-UC-PI-Y-P	GRA5234
Pre-Inflated	Youth	10	PHPB-UC-PI-Y-P5	GRA5235
Pre-Inflated	Small	1	PHPB-PI-SM-P	GRA5228
Pre-Inflated	Small	10	PHPB-PI-SM-P5	GRA5229
Pre-Inflated	Standard	1	PHPB-UC-PI-ST-P	GRA5230
Pre-Inflated	Standard	10	PHPB-UC-PI-ST-P5	GRA5231
Pre-Inflated	X-Large	1	PHPB-UC-PI-XL-P	GRA5232
Pre-Inflated	X-Large	10	PHPB-UC-PI-XL-P5	GRA5233

Sizing Guide

Size**	Shoe Size UK	Shoe Size EU
Child	c5 - c8	22 - 26
Youth	c9 - 1	27 - 33
Small	2 - 7	34 - 40
Standard	7.5 - 11	41 - 45
X-Large	11.5+	46+

^{**} If calf circumference is larger than normal (maybe oedema present) consider fitting larger size than determined by shoe size.





HeelSafe

Pressure Redistribution Pad

- Pre-inflated
- A solution for ambulatory patients and those at risk of falls when wearing boots
- An alternative choice for patient refusing or unsuitable for boots
- Multi-patient use
- Pressure mapped⁴
- Can be combined with SoleSafe
- One size 80 x 66cm (W x L)
- * Supplied individually
- * Multi patient use


Ideal for patients at risk of a fall

Style	Size	Pump	Size WxL (cm)	Product Code
Covered	Pre-Inflated		80 x 66	PHP-C-PI-OS
Covered	Inflatable	Yes	80 x 66	PHP-C-IP-OS
Covered	Inflatable	No	80 x 66	PHP-C-NP-OS

SoleSafe

Bed End Pressure Redistribution Pad

- Pre-inflated
- Protects the soles of feet
- Allows objects hung on the end of the bed to be concealed
- Pressure mapped³
- Can be combined with HeelSafe above
- One size 80 x 28cm (W x L)
- * Supplied individually
- * Multi patient use





Ideal for patients who slide down the bed

Footshield II

Shoe

- Evidence based rocker sole⁵
- 10° toe out angle
- Apex of rocker at 60% of shoe length
- 20° rocker angle
- Removable memory foam insole
- Square toe design accommodates K-Wires
- Straps over bulky dressings
- Can be used with optional Hex offloading insole on page 8
- Removable toe cap



-1-				•		
~	Supp	חסוו	Ind	111/10	liial	1
	Jupp	II E U	HIIU	IIVIU	ıuaı	ıν
						•

" Supplied individually								
				LEFT			RIGHT	
Size	UK Shoe Size	Length (mm)	Product Code	NHSSC Code	PIP Code	Product Code	NHSSC Code	PIP Code
X-Small	1.5 - 3	240	FSS-II/XSL	ERP5998		FSS-II/XSR	ERP5999	
Small	3.5 - 5	257	FSS-II/SL	ERP5994	422-5991	FSS-II/SR	ERP5995	422-6296
Medium	5.5 - 7	273	FSS-II/ML	ERP5992	422-6478	FSS-II/MR	ERP5993	422-6502
Large	7.5 - 9	290	FSS-II/LL	ERP5990	422-6270	FSS-II/LR	ERP5991	422-6288
X-Large	9 - 11	307	FSS-II/XLL	ERP5996	422-6486	FSS-II/XLR	ERP5997	422-6494
WITH HEX	(INSOLE							
X-Small	1.5 - 3	240	FSS-IIH/XSL	ERP5988		FSS-II/XSR	ERP5989	
Small	3.5 - 5	257	FSS-IIH/SL	ERP5984		FSS-IIH/SR	ERP5985	
Medium	5.5 - 7	273	FSS-IIH/ML	ERP5982		FSS-IIH/MR	ERP5983	
Large	7.5 - 9	290	FSS-IIH/LL	ERP5980		FSS-IIH/LR	ERP5981	
X-Large	9 - 11	307	FSS-IIHH/XLL	ERP5986		FSS-IIH/XLR	ERP5987	

Footshield II

Boot

• Features as above for the Footshield II Shoe

* Supplied individually

		LEFT		RIGHT			
Size	UK Shoe Size	Length (mm)	Product Code	NHSSC Code	Product Code	NHSSC Code	
X-Small	1.5 - 3	240	FSB-II/XSL	ERP6018	FSB-II/XSR	ERP6019	
Small	3.5 - 5	257	FSB-II/SL	ERP6014	FSB-II/SR	ERP6015	
Medium	5.5 - 7	273	FSB-II/ML	ERP6012	FSB-II/MR	ERP6013	
Large	7.5 - 9	290	FSB-II/LL	ERP6010	FSB-II/LR	ERP6011	
X-Large	9 - 11	307	FSB-II/XLL	ERP6016	FSB-II/XLR	ERP6017	
WITH HEX INSOLE							
X-Small	1.5 - 3	240	FSB-IIH/XSL	ERP5978	FSB-II/XSR	ERP5979	
Small	3.5 - 5	257	FSB-IIH/SL	ERP6034	FSB-IIH/SR	ERP5975	
Medium	5.5 - 7	273	FSB-IIH/ML	ERP6032	FSB-IIH/MR	ERP6033	
Large	7.5 - 9	290	FSB-IIH/LL	ERP6030	FSB-IIH/LR	ERP6031	
X-Large	9 - 11	307	FSB-IIH/XLL	ERP5976	FSB-IIH/XLR	ERP5977	



Square toe design to accommodate K-Wires





Footshield II

Wedge Shoe

- Forefoot off-loader
- Removable memory foam insole
- Square toe design accommodates K-Wires
- Straps over bulky dressings
- Can be used with optional Hex offloading insole below
- Removable toe cap
- * Supplied individually

			LEFT		RIGHT		
Size	UK Shoe Size	Length (mm)	Product Code	NHSSC Code	Product Code	NHSSC Code	
X-Small	1.5 - 3	240	FSWS-II/XSL	ERP6028	FSWS-II/XSR	ERP6029	
Small	3.5 - 5	257	FSWS-II/SL	ERP6024	FSWS-II/SR	ERP6025	
Medium	5.5 - 7	273	FSWS-II/ML	ERP6022	FSWS-II/MR	ERP6023	
Large	7.5 - 9	290	FSWS-II/LL	ERP6020	FSWS-II/LR	ERP6021	
X-Large	9 - 11	307	FSWS-II/XLL	ERP6026	FSWS-II/XLR	ERP6027	
WITH HEX INSOLE							
X-Small	1.5 - 3	240	FSWS-IIH/XSL	ERP5988	FSWS-IIH/XSR	ERP5989	
Small	3.5 - 5	257	FSWS-IIH/SL	ERP5984	FSWS-IIH/SR	ERP5985	
Medium	5.5 - 7	273	FSWS-IIH/ML	ERP5982	FSWS-IIH/MR	ERP5983	
Large	7.5 - 9	290	FSWS-IIH/LL	ERP5980	FSWS-IIH/LR	ERP5981	
X-Large	9 - 11	307	FSWS-IIH/XLL	ERP5986	FSWS-IIH/XLR	ERP5987	

Removeable memory foam insole

Footshield

Hex Insole

- Allows for off-loading of localised wounds and ulcers
- Adhesive stabiliser board maintains integrity of void
- Removable hexagonal pegs create a void
- * Supplied individually



			SHOE/BOOT		WEDGE SHOE	
Size	UK Shoe Size	Length (mm)	Product Code	NHSSC Code	Product Code	NHSSC Code
X-Small	1.5 - 3	240	FSHISB/XS	ERF616	FSHIWS/XS	ERF614
Small	3.5 - 5	257	FSHISB/S	ERF617	FSHIWS/S	ERF612
Medium	5.5 - 7	273	FSHISB/M	ERF618	FSHIWS/M	ERF611
Large	7.5 - 9	290	FSHISB/L	ERF619	FSHIWS/L	ERF610
X-Large	9 - 11	307	FSHISB/XL	ERF615	FSHIWS/XL	ERF613

Ambulatory Pressure

Relief Boot

- Offloads pressure from heel
- Soft, breathable and washable liner
- Adjustable length toe support
- Anti-rotation bar
- Available in two sizes

* Supplied individually

Size	Calf Circ. (cm)	Order Code	NHSSC Code
Standard	Max. 45	APRB/ST	GRD3625
Large	38 - 55	APRB/L	GRD3624



Pressure Relief

AFO

- Universal size
- Offloads pressure from heel
- Tri-laminate material allows air to circulate and helps wick moisture
- Adjustable length toe support
- Padded ankle strap reinforces foot position
- Supplied with spare liner set
- Non-slip sole allows for ambulation
- * Supplied individually

Size	Product Code	NHSSC Code
Universal	AUPRNS/U	GRD3627





AT RISK PATIENTS

Product choice pathway

AMBULANT PATIENT

Fit device to bed





NON-AMBULANT PATIENT

Patient worn or fit to bed









ACTIVE FOOT OR HEEL ULCER PATIENTS Product choice pathway AMBULANT PATIENT NON-AMBULANT PATIENT

Rigid frame device





Non-rigid device





How to Order:

Ordering any product within our pressure care collection is easy. You can contact us by phone, email or online using our contact details below.

If you need help selecting products or have any questions, please get in touch, our Customer Services Team will be happy to assist you.



01246 268456



info@talarmade.com



talarmade.com

PRESSURE CARE PRODUCT RANGE


Products designed for use in the in the prevention and management of pressure related injuries.

















Made with clinical excellence in mind

Edinburgh House Millennium Way Dunston Technology Park Chesterfield S41 8ND



01246 268456



info@talarmade.com



talarmade.com

References

- Reducing PU incident rates by giving the run around. Hayley Jones, Tissue Viability Advanced Nurse Specialist & Gemma McGrath, Tissue Viability Clinical Nurse Specialist, Medway NHS Foundation Trust. Poster, Wounds UK Annual Conference, November 2017.
- Evaluation of a pressure ulcer prevention offloading boot in the acute setting. Deborah Gleeson, Lead Nurse, Tissue Viability. August 2017.
- 3. Evaluation of Prolevo SoleSafe. Duncan Bain Consulting, available on request.
- 4. Evaluation of Prolevo HeelSafe. Duncan Bain Consulting, available on request.
- Effect of rocker sole design features on forefoot plantar pressures in people with and without diabetes. Chapman et al, Clinical Biomechanics, 28(6) 679-85, July 2013.