

EGN08 - NITREX Extra Length

Blue nitrile powder free examination gloves with a long cuff



Latex Free

Powder Free

Micro Textured Fingers

Extended Cuff for Additional Coverage

Specification

Material	Nitrile
Length	300mm ± 10mm
Protein Level	Contains no latex proteins
Surface	Micro Textured Fingers
Shape	Ambidextrous
Colour	Blue
Sterilisation	Non sterile
Shelf Life	3 years from manufacture
Packing	1 case = 1000 gloves (10 dispensers of 100)

Quality standards: Manufactured in a facility holding ISO 9001 and ISO 13485.

Product standards:

Medical Device

This product is classified as a Class I Medical Device and is tested to European Standards EN 455-1:2020, EN 455-2:2015, EN 455-3:2015 and EN 455-4:2009.

<u>PPE</u>

This product is classed as Category III Personal Protective Equipment (PPE) according to the PPE Regulation (EU) 2016/425. EU Type Examination carried out by SATRA. Tested to European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 and EN ISO 374-5:2016.

Re-Order Codes

Extra Small	EGN08E
Small	EGN08S
Medium	EGN08M
Large	EGN08L
Extra Large	EGN08X



Storage: Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

Country of origin: Malaysia

Notice: This product is free from natural rubber latex and is suitable for use in latex-free environments. Contains accelerators, which may cause allergic reactions. Stop use of the product if skin reactions occur. Single use only.

Find out more at www.medicareproducts.com