







# **BioTouch (GM008\*)**

biodegradable · AQL 1.5 · chemical tested · violet blue medical grade · medium weight · nitrile · examination gloves



# **Physical Specifications:**

Material: Nitrile (latex free) Colour: Violet blue Length: 240 mm Weight: 3.5 g

Thickness: 0.07 mm Finish: Powder free Fit: Ambidextrous Sterile: No

## **Order Codes and Size** Information:

GM0082 | Small GM0083 | Medium GM0084 | Large GM0085 | Extra large

#### **Packaging:**

Each dispenser box contains 100 gloves Each case contains 10 dispenser boxes

#### Classification:









For more about the Advanced range, please click on the



# **⇔** Features

- 90% biodegradation in 490 days\*
- Contains innovative biodegradable technology
- Tested to EN374
- Tested to ASTM D5511 and ASTM D5526
- AQL 1.5 and EN455 parts 1-4
- Tested to EN1186
- Nitrile formulation conforms to the hand during use
- Manufactured with textured fingertips
- Touchscreen compatible

### **⇔** Benefits

- Environmentally friendly, much faster organic break down compared to regular nitrile, which can take over 100 years
- The technology means the biodegradation process only starts when disposed of in an active landfill environment
- Resistant to a range of chemicals
- Tested to leading biodegradation standards
- Meets medical standards including pinholes
- Suitable for contact with all food types
- Nitrile formulation conforms to the hand during use for a more comfortable fit
- Beaded cuff provides extra strength during donning and helps prevent liquid roll back
- Textured fingertips provide improved grip, particularly when working in wet or oily conditions
- Suitable for use with mobile and tablet displays without removing gloves

## **⇔** Sectors

- Food
- Automotive
- Aerospace
- Chemical manufacture and use
- Dentistry
- Dental hygiene
- Electronics
- Engineering
- Food processing
- Healthcare
- Laboratory work
- · Manufacturing and assembly
- Pharmaceutical work
- Printing



\* Above results are based on a real-time study. The actual biodegradation rates will vary depending on the landfill conditions and the biological activity of microorganisms surround the gloves.

#### UNIGLOVES UK

Registered in England No: 4010200 VAT registration 791817787

- 3 Ambley Green, Gillingham Business Park, Gillingham, ME8 ONJ, UK
- +44 (0)800 049 6602
- enquiries@unigloves.co.uk unialoves.co.uk

# UNIGLOVES GERMANY

Arzt- & Klinikbedarf Handelsgesellschaft mbH

- Camp-Spich-Straße 71 D-53842 Troisdorf-Spich,
- +49 (0) 2241-9323-0 unigloves@unigloves.de unialoves.de



fingertip

Issue No. 4

Prepared at Gillingham on 08/03/22