

## GN01 - NITREX Extra Sensitive

Blue nitrile powder free examination gloves.



Latex Free

Powder Free

Micro Textured Fingers

Ultimate Tactile Sensitivity

### Specification

<b>Material</b>	Nitrile
<b>Length</b>	240mm
<b>Protein Level</b>	Contains no latex proteins
<b>Surface</b>	Micro Textured Fingers
<b>Shape</b>	Ambidextrous
<b>Colour</b>	Blue
<b>Sterilisation</b>	Non sterile
<b>Shelf Life</b>	3 years from manufacture
<b>Packing</b>	1 case = 2000 gloves (10 dispensers of 200)

**Quality standards:** Manufactured in a facility holding ISO 9001 and ISO 13485.

#### Product standards:

##### Medical Device

This product is classified as a Class I Medical Device and is tested to European Standards EN 455-1:2020, EN 455-2:2015, EN 455-3:2015 and EN 455-4:2009.

##### PPE

This product is classed as Category III Personal Protective Equipment (PPE) according to the PPE Regulation (EU) 2016/425.

EU Type Examination carried out by SATRA. Tested to European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016 and EN ISO 374-5:2016.

### Re-Order Codes

Extra Small	GN01E
Small	GN01S
Medium	GN01M
Large	GN01L
Extra Large	GN01X



**Storage:** Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

**Country of origin:** Malaysia

**Notice:** This product is free from natural rubber latex and is suitable for use in latex-free environments. Contains accelerators, which may cause allergic reactions. Stop use of the product if skin reactions occur. Single use only.

Find out more at [www.medicareproducts.com](http://www.medicareproducts.com)