

GN06 - NITREX Accelerator Free

Indigo blue nitrile examination gloves.



Latex Free

Powder Free

Accelerator Free

Enhanced Tactile Sensitivity

Specification

Material	Nitrile
Length	240mm
Protein Level	Contains no latex proteins
Surface	Micro Textured Fingers
Shape	Ambidextrous
Colour	Indigo blue
Sterilisation	Non sterile
Shelf Life	3 years from manufacture
Packing	1 case = 2000 gloves (10 dispensers of 200)

Quality standards: Manufactured in a facility holding ISO 9001 and ISO 13485.

Product standards:

Medical Device

This product is classified as a Class I Medical Device and is tested to European Standards EN 455-1:2020, EN 455-2:2015, EN 455-3:2015 and EN 455-4:2009.

PPE

This product is classed as Category III Personal Protective Equipment (PPE) according to the PPE Regulation (EU) 2016/425. EU Type Examination carried out by SATRA. Tested to European Standards EN 420:2003+A1:2009, EN ISO 374-1:2016+A1:2018 and EN ISO 374-5:2016.

Re-Order Codes

Extra Small	GN06E
Small	GN06S
Medium	GN06M
Large	GN06L
Extra Large	GN06X



Storage: Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

Country of origin: Malaysia

Notice: This product is free from natural rubber latex and is suitable for use in latex-free environments. Stop use of the product if skin reactions occur. Single use only.

Find out more at www.medicareproducts.com