

# **GV05 - Medisafe Vytrex**

Clear vinyl powder free examination gloves.



**Latex Free** 

**Powder Free** 

**Easy Donning** 

**Smooth Surface** 

## **Specification**

**PVC Plastisol Material** 

Length 240mm

**Protein Level** Contains no latex proteins

**Surface** Smooth

Ambidextrous Shape

Colour Clear

Sterilisation Non sterile

**Shelf Life** 5 years from manufacture

**Packing** 1 case = 1000 gloves (10 dispensers of 100)

Quality standards: Manufactured in a facility holding ISO 9001 and ISO 13485.

#### **Product standards:**

#### **Medical Device**

This product is classified as a Class I Medical Device and is tested to European Standards EN 455-1:2020, EN 455-2:2015, EN 455-3:2015 and EN 455-4:2009.

### **Re-Order Codes**

Small **GV05S** GV05M Medium GV05L Large GV05X Extra Large















Storage: Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

Country of origin: China

**Notice:** This product is free from natural rubber latex and is suitable for use in latex-free environments. Stop use of the product if skin reactions occur. Contains accelerators which may cause allergic reactions. Single use only.