


## Overchair Table - Height Adjustable, Melamine

### Product Images



### Description

#### FEATURES

- BS 2483: 1977 Specification for overbed tables (selected tests)
- BS EN 438-2: 2016 High-pressure decorative laminates (HPL). Sheets based on thermosetting resins (usually called laminates). Determination of properties
- Height adjustable using a gas damper (735 - 1045mm)
- Safety feature prevents trapping of the patient between the bed / chair & table
- Extruded aluminium column provides a smooth linear height adjustment
- Cantilever tubular base available as...
  - C shaped suitable for use with beds & chairs
- Melamine top with varnished edges available with either a flat top or one varnished beech lip
- 50mm polypropylene / thermoplastic

#### OVERALL DIMENSIONS (W X D X H)

- Overchair Table - 918 x 440 x 735 - 1045mm
- Top - 880 x 395mm

STANDARD COLOURS

- Top - Grey White / Beech
- Upper Extrusion - Spearmint / Grey White
- Lower Extrusion / Base - Grey White

Code	Description	Overall Dimensions (W x D x H)	Weight
OCL/GW0/WH/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/GW/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Grey White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/SP/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Spearmint Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/FR/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Flame Red Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/JB/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Jet Black Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/SY/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Signal Yellow Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/LB/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Light Blue Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/OR/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Orange Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/PU/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Purple Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW0/SG/S	Overchair Table - Height Adjustable, Melamine, Grey White Flat Top with Signal Green Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/WH/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/GW/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Grey White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/SP/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Spearmint Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/FR/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Flame Red Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/JB/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Jet Black Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/SY/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Signal Yellow Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/LB/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Light Blue Upper Extrusion	918 x 440 x 735-1045mm	14.1kg

OCL/BE0/OR/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Orange Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/BE0/PU/S	Overchair Table - Height Adjustable, Melamine, Beech Flat Top with Purple Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW1/WH/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW1/GW/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Grey White Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW1/SP/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Spearmint Upper Extrusion	918 x 440 x 735-1045mm	14.1kg
OCL/GW1/FR/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Flame Red Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/JB/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Jet Black Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/SY/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Signal Yellow Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/LB/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Light Blue Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/OR/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Orange Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/PU/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Purple Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/GW1/SG/S	Overchair Table - Height Adjustable, Melamine, Grey White Top with One Lip, with Signal Green Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/WH/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with White Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/GW/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Grey White Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/SP/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Spearmint Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/FR/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Flame Red Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/JB/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Jet Black Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/SY/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Signal Yellow Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/LB/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Light Blue Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/OR/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Orange Upper Extrusion	918 x 440 x 735-1045mm	15kg
OCL/BE1/PU/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Purple Upper Extrusion	918 x 440 x 735-1045mm	15kg

OCL/BE1/SG/S	Overchair Table - Height Adjustable, Melamine, Beech Top with One Lip, with Signal Green Upper Extrusion	918 x 440 x 735-1045mm	15kg
--------------	--	------------------------	------

### Additional Information

---

Colour	Brown
--------	-------