

FALL SAVERS

Nurse Call System List

VARIOUS CONNECTORS

¼" Mono Jack



¼" Stereo Jack



3.5mm Jack



4P4C RJ (RJ 9)



6P6C RJ (RJ11) Rev BT















3 Pin Din








BNC













Nurse Call System	Connector Type	Stock Adapter Available (Y/N)	Adapter Part #
<p>Aidcall ATX 3000</p> 	<p>Reversed BT plug</p> 	<p>Yes</p>	<p>Adapter-BT</p>
<p>Aidcall ATX 3000+/4000/5000</p> 	<p>6P6C RJ</p> 	<p>Yes</p>	<p>Adapter-ATX6</p>
<p>Attendo</p>	<p>¼" mono jack (1 ring on jack plug)</p> 	<p>N/A</p>	<p>N/A</p>

<p>ARM Nursecall 2000</p> 	<p>¼" stereo jack (2 rings on Jack plug)</p> 	<p>N/A</p>	<p>N/A</p>
<p>Blick</p>	<p>¼" stereo jack (2 rings on Jack plug)</p> 	<p>N/A</p>	<p>N/A</p>
<p>C-Tec Nursecall (600/800 series)</p> 	<p>¼" stereo jack (2 rings on Jack plug)</p> 	<p>N/A</p>	<p>N/A</p>
<p>Chubb</p> 	<p>¼" mono jack (1 ring on jack plug)</p> 	<p>N/A</p>	<p>N/A</p>

<p>Courtney Thorne</p> 	<p>¼" mono jack with diode (not standard ¼" mono jack)</p> 	<p>Yes</p>	<p>Adapter-CT</p>
<p>Eclipse</p> 	<p>¼" stereo jack (Tip & Ring connection)</p>	<p>No</p>	<p><i>"Jack socket/plug to ¼" stereo jack (Tip & Ring connection)"</i></p> <p>Custom Lead</p>
<p>Intercall</p> 	<p>¼" mono jack (1 ring on jack plug)</p> 	<p>N/A</p>	<p>N/A</p>

<p>LiveLink (ECO1-11/ECO1-12)</p> 	<p>4P4C RJ</p> 	<p>No</p>	<p><i>"Jack socket/Plug to 4pin RJ (wired to pins 1 & 2)"</i></p> <p>Custom Lead</p>
<p>Mawtec (Air-Aid scan call)</p>	<p>3 Pin Din</p> 	<p>No</p>	<p><i>"Jack socket/Plug to 3 pin din"</i></p> <p>Custom Lead</p>
<p>Medicare (1000 & 4000 series)</p> 	<p>4P4C RJ</p> 	<p>Yes</p>	<p>Adapter-MED4</p>

<p>Medicare (6000 series)</p> 	<p>6P6C RJ</p> 	<p>Yes</p>	<p>Adapter-MED6</p>
<p>Quantec</p> 	<p>¼" stereo jack (2 rings on Jack plug)</p> 	<p>N/A</p>	<p>N/A</p>
<p>Scantronic</p>	<p>BNC plug</p> 	<p>No</p>	<p><i>"Jack socket/Plug to BNC plug"</i></p> <p>Custom Lead</p>
<p>SAS Network 2</p> 	<p>¼" stereo jack (body/ring connection) – please enquire</p> 	<p>No</p>	<p>Please enquire</p> <p>Custom Lead</p>

<p>Targon Technologies</p>	<p>¼" mono jack (1 ring on jack plug)</p> 	<p>N/A</p>	<p>N/A</p>
<p>Tynetec</p> 	<p>3.5mm plug (screw type)</p> 	<p>No</p>	<p><i>"Jack socket/Plug to 3.5mm plug"</i></p> <p>Custom Lead</p>